



AN EFFECTIVE MEDIUM AGAINST DERMATOLOGICAL DISORDERS:



Waldmann is one of the world's leading manufacturer of diagnostic and therapy systems for dermatological phototherapy. Based on decades of long experience and findings gained from the continuous, intensive dialogue with dermatologists, the criteria for forward-thinking system concepts were analyzed and innovatively implemented.

As a whole, but also with respect to the multifaceted details, the result marks the technical progress, which creates new possibilities, offers interesting economical perspectives, and is compatible with state-of-the-art interactive systems.

Safety has first priority, and in addition to the expanded treatment possibilities, new technical equipment must prominently fulfill this safety aspect.


A key aspect is the realization that illnesses create uncertainty and generate fear. As a result, disease patterns that can be indicated by UV therapies require unconditional mutual trust between the treating physician and the patient.

Waldmann therapy systems promote this process through diverse, exactly meterable therapy protocols, and wellengineered safety technology.

This is something you can rely on.

The treatment success, the associated intensity of patient commitment, as well as the reduction of economic risks are closely related to safety and trust. The team of Waldmann experts supports the dermatologist in all areas and optimizes the prospects for success.

UV TECHNOLOGY FROM WALDMANN



UNIFORM FULL BODY UV THERAPY

The professional full body therapy systems from Waldmann can be used for monotherapy and combination therapy. All appliances with dosage entry in Joules/cm².

UV 7002

Convenient to operate therapy system with constant measurement and dosage adjustment. Extra long UV lamps. Maximum safety due to control via the practice's PC with therapy software.

UV 5040

Compact full body UV therapy system. Can be fitted with various lamps.

UV 1000

Total body irradiation cubicle with LCD timer, automatic conversion of the Joule value (J/cm²) into treatment duration. Low space requirements.

UV 100 L

For use as a monotherapy appliance in a physician's office or for use in home therapy. Can be fitted with various lamps. Fold-out sides for irradiation of the sides of the body. Low space requirements.



UV 7002



UV 5040



UV 1000



UV 100 L



PARTIAL THERAPY SYSTEMS: MOBILE, VERSATILE AND EFFICIENT



VARIOCONTROL

VARIOCONTROL
Important measurement accessory for determining the irradiation intensity (mW/cm²) for systems without built-in sensors. A complete skin test system when used in conjunction with the TEST UNIT (optional accessory).



TEST UNIT

TEST UNIT
Optional accessory to the VARICON-TROL device for the determination of erythema doses, separated according to UVA and UVB.



UV 109
WITH ADD-ON COMB

UV 109
Hand-held therapy system with UVB broad band, UVB narrow band or UVA lamp. With add-on comb for separating hair.

UV 182
Flexible system for a variety of partial treatments with UVA, narrowband UVB or broadband UVB lamps. Modules with different lamps can be easily inserted into the basic device and swiveled to the required position. With pivoting TFT-Monitor with resistive touch screen for central control of the therapy modules. Automatic module recognition.



UV 182



UVA 700 L

UVA 700 L
High performance system for photoprovocation tests and therapeutic UVA and PUVA applications. High dosage ranges possible.

UV 802 L
Partial body therapy systems for versatile use. Monotherapy: UVA, PUVA, UVB narrow band, UVB broad band. Combination therapy: UVA and UVB narrow band synchronously, UVA and UVB broad band synchronously. For diagnosis: MPD and MED tests with VARIOCONTROL and TEST UNIT. With LCD timer for dosage entry in Joules/cm².



UV 802 L



DIAGNOSIS SYSTEMS

A reliable diagnosis is the base for every successful therapy. Waldmann offers a full variety of diagnosis systems in different versions and with many optional accessories.



DHL 111 M

DHL 111 M
Easy to handle and lightweight. With compact fluorescent tube.

LRE 122 T
Diagnosis light on stand. Multiple mounting options: Wall bracket, table foot or clamp, 5-feet on castors.



LRE 122 T



DL 122 M

DL 122 M
Diagnosis light for versatile use in all dermatological examinations.

DHL 109 M
Handy and lightweight hand wood-light for detection of fluorescence of disorders and assessment of pigment changes.



DHL 109 M

DHLL 404 M / DHLL 204 UM
Single Wood's light or as combination model with a toggle switch for neutral white fluorescent tube light.



DHLL 404 M / DHLL 204 UM

Wall bracket available as accessory.

PHOTODYNAMIC THERAPY

Waldmann, leading manufacturer with more than 30 years of experience in phototherapy decisively influenced the development of Photodynamic Diagnosis and Therapy in cooperation with the most important clinical centres in Europe.

PDT 1200 L

The PDT 1200 L illuminates an area of up to 15 cm diameter homogeneously and with a high irradiance.

The freely adjustable irradiance of up to 200 mW/cm² may especially be important for the treatment of oncologic indications.



PDT 1200 L



CUSTOMER SERVICE AND MAINTENANCE

Waldmann medical systems are developed and produced according to the most up-to-date technology.

As a leading manufacturer of therapy systems, Waldmann also guarantees a first-rate customer service.

In many cases, faults can be remedied simply by troubleshooting. Waldmann offers favourable rates and conditions for onsite customer service. A maintenance contract includes regular checks and guarantees smooth running of your therapy system.

Herbert Waldmann GmbH & Co. KG
 Peter-Henlein-Straße 5
 D-78056 VILLINGEN-SCHWENNINGEN

FREE INFORMATION BY POST OR FAX

PLEASE SEND ME AN OFFER FOR:

Model

Lamp configuration

INFORMATION ON THE TOPIC:

- | | |
|---|--|
| <input type="checkbox"/> Professional UV therapy systems | <input type="checkbox"/> Diagnosis Systems |
| <input type="checkbox"/> Waldmann UV home Therapy Systems | <input type="checkbox"/> Lamps |
| <input type="checkbox"/> Photodynamic therapy | <input type="checkbox"/> Other _____ |

- I would like to make an appointment with your medical product consultant.
 Please contact:

Phone _____

Clinic

Mr./Ms./Prof./Dr.

Department

Street

Post code, location

e-Mail

FAX NO. : + 49 7720 601-425

HEADQUARTERS GERMANY

Herbert Waldmann GmbH & Co. KG
Postfach 50 62
78057 VILLINGEN-SCHWENNINGEN
GERMANY
Telephone +49 7720 601-200
Fax +49 7720 601-290
www.waldmann-medizintechnik.com
info-puva@waldmann.com

FRANCE

Waldmann Eclairage S.A.S.
Zone Industrielle
Rue de l'Embranchement
67116 REICHESTETT
FRANCE
Telephone +33 3 88 20 95 88
Fax +33 3 88 20 95 68
www.waldmann.com
info-fr@waldmann.com

SWITZERLAND

Waldmann Lichttechnik GmbH
Benkenstrasse 57
5024 KÜTTIGEN
SWITZERLAND
Telephone +41 62 839 1212
Fax +41 62 839 1299
www.waldmann.com
info-ch@waldmann.com

THE NETHERLANDS

Waldmann BV
Lingewei 19
4004 LK TIEL
THE NETHERLANDS
Telephone +31 344 631019
Fax +31 344 627856
www.waldmann.com
info-nl@waldmann.com